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
Component

English

※ Optional replacement heads (two-row, angled, regular) sold separately

Handle

Wireless charging cradle

Replacement brush head storage case

User manual

1

Product structure

English

Toothbrush head

Power/Action button

Mode indicator LED


Charging indicator LED

Wireless charging cradle


2

Product function

English

Replacing the Brush Head

- **Installation:** Install the brush head as shown in the illustration on the left. Ensure the bristles are facing forward. Do not use excessive force if it is not aligned correctly, as this may damage the product.
- **Included Items:** Replacement brush heads are included in the package.
- **Replacement Interval:** We recommend replacing the brush head every 3 months. Replace it sooner if the bristles become worn or deformed. (Additional brush heads are sold separately.)
- **Compatibility:** This device is only compatible with genuine brush heads designed for this model. Using incompatible third-party heads may cause the device to malfunction.
- **1mm Gap:** A gap of approximately 1mm between the handle and the brush head is normal and required for proper vibration. Do not attempt to close this gap by pressing too hard, as it may cause damage to the device.


3

How to use

English

- 1. Preparation**
Apply an appropriate amount of toothpaste to the bristles. To prevent splashing, the intensity increases gradually over the first 3 seconds. Start brushing within this 3-second window.
- 2. Power ON & Intensity Control**
Press the power button briefly to turn the device on. While on, press the button again briefly to cycle through intensity levels. Select the intensity that suits you best.
- 3. Brushing & Power OFF**
A 2-minute brushing session is recommended. The device will automatically turn off after 2 minutes. To turn off the power manually during use, press and hold the power button for more than 1 second.
- 4. Cleaning & Storage**
After use, rinse the bristles thoroughly under running water and let them air dry. Store the device in a cool, dry place.

Smart 2-Minute Timer · 30-Second Interval Alert

- This product features a built-in timer to guide your brushing time.
- The session lasts for 2 minutes, after which the power shuts off automatically.
- To ensure even cleaning, the mouth is divided into 4 quadrants. Every 30 seconds, the device will pause briefly to signal you. Move to the next quadrant whenever you feel this vibration/pause.

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Mode Settings

- bright Implant Care features three modes: Gum, Daily, and Intensive. Choose the mode that best suits your needs.
- To switch modes, briefly press the power button while the device is ON.The modes will cycle in the following order: Gum → Daily → Intensive.
- Memory Function:The device remembers your last used mode and will automatically start with that mode during your next session.

Power/Mode Button

Intensive

Daily

Gum

Cleaning & Storage

- Rinse the brush head thoroughly under running water immediately after use and allow it to dry.
- Residual toothpaste or improper drying can lead to bacterial growth.
- Once cleaned and dried, store the brush head in a cool, dry place.

- Lightly clean the handle with water or a damp cloth after use.
- Wipe with a dry towel and allow it to air dry.
- As this product contains a built-in battery, do not store it in hot or humid environments.

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Tel. 031-888-5100 Fax. 031-8019-9123

Charging

- **Battery Life:** Lasts up to 50 days on a full charge.(Based on brushing twice daily for 2 minutes: 30 days on High, 40 days on Medium, 50 days on Low)Battery life may vary depending on usage.
- **Wireless Charging:** Supports wireless charging. Place the handle on the powered charging base to begin charging.
- **Charging Indicator:** Red light (below 80%), blue light (80% or above).
- **Usage Indicator:** Red light (below 20%), blue light (20% or above).
- **Safety Feature:** The device will not operate when placed on or near the charging base.
- **Power Specifications:** Rated at 5V 1A. Use a compatible adapter.(Charging may not work with adapters below 5V 1A.)
- **Warning:** Ensure the USB power connection is dry. Moisture may cause fire or electric shock.

Maintenance & Support

Warranty (A/S)

This product includes a 6-month limited warranty. After 6 months, repair services will be subject to a fee.

For service inquiries, please contact the seller or manufacturer.

Note: Damage caused by user misuse or negligence is not covered under the warranty.

Disposal

This product contains a built-in battery and electronic components. Do not dispose of with household waste. Please dispose of it at a designated electronic waste collection point to help protect the environment.

Precautions

User Safety

- Read all instructions carefully before use.
- Consult a dentist if you have recently undergone gum surgery or have severe gum swelling.
- Stop use and consult a doctor if you experience severe pain or discomfort.
- Keep out of reach of children under 3 years of age. Children under 13 or individuals requiring assistance should use this product under supervision.
- Do not apply excessive pressure, as it may damage teeth or gums or cause product malfunction.

Charging Safety

- Charge the device when the red LED indicator at the bottom flashes.
- This product contains a built-in battery. Avoid exposure to high temperatures or strong impacts to reduce the risk of fire or explosion.
- Use a compatible 5V/1A (5W) adapter. Using incompatible chargers may cause malfunction.
- Do not use the charger near water (e.g., bathtubs or sinks). Ensure the USB port is completely dry before charging to prevent fire or electric shock.

Maintenance Precautions

- Use only genuine replacement brush heads. Other types are not compatible.
- Rinse the brush head thoroughly after use and allow it to air dry to prevent bacterial growth.
- Replace the brush head every 3 months or sooner if the bristles become deformed.
- Do not immerse the handle or brush head in boiling water, detergents, or chemicals.
- Store the product in a cool, dry place. Exposure to extreme heat or humidity may cause damage.

Warranty (A/S)

This product includes a 6-month limited warranty. After 6 months, repair services will be subject to a fee.

For service inquiries, please contact the seller or manufacturer.

Note: Damage caused by user misuse or negligence is not covered under the warranty.

Disposal

This product contains a built-in battery and electronic components. Do not dispose of with household waste. Please dispose of it at a designated electronic waste collection point to help protect the environment.

FAQ

Q. The vibration is too weak.

→ Applying excessive pressure can reduce vibration. Use gentle pressure.

→ Attaching the brush head too tightly may affect motor performance. Ensure there is a 1 mm gap between the head and handle. If there is no gap and the vibration remains weak, contact Customer Service.

Q. It doesn't work even after charging.

→ The battery may be completely discharged, or the power button may not have been pressed correctly. Fully charge the device and try again.

→ The battery may have reached the end of its lifespan. If the device does not charge or operate after a full charge, contact Customer Service.

Q. The vibration is too strong / My gums hurt after use.

→ If you have sensitive gums, the vibration may feel strong. Use a lower intensity mode, or discontinue use and consult a dentist.

Q. It gets very hot during charging.

→ Wireless charging may generate some heat. However, if the device becomes excessively hot, disconnect the charger immediately and contact Customer Service.

Q. It makes a strange noise during operation.

→ Some noise is normal for motorized devices. If you hear continuous or unusual noise, contact Customer Service.

Q. Can I leave the fully charged handle on the charger?

→ The device includes overcharge protection. However, to extend battery life, we recommend removing it from the charger once fully charged.

Q. Can I keep the charger and handle in a humid bathroom?

→ The handle is rated IPX8 and is waterproof. However, the charger is not waterproof. Store both in a cool, dry place to prevent moisture buildup.

Q. Can I wash the handle under running water?

→ Yes, the handle is rated IPX8 and can be rinsed under running water. However, do not wash the charger, as it is not waterproof.

Product Specifications

Product Name: bright Implant Care

Models: BSET-001, BSET-002, BSET-003, BSET-004, BSET-005, BSET-006, SIP-WC-1.5

Materials: ABS, PC, POM, etc.

Rated Input: DC 5V, 1A

Power Consumption: 2 W

Motor Vibration: 21,000 vpm (vibrations per minute)

Wireless Charging Frequency: 100 kHz ± 10%

Battery Specifications: 3.7 V, 1100 mAh

Waterproof Rating: IPX8

Dimensions: 180 × 28.5 × 27.5 mm (handle)

Weight: 116 g

Manufacturer information

Manufacturer: GENOSS Co., Ltd. (South Korea)

D-factory, 56, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, Republic of Korea

Tel: +82-31-888-5100 | Web: www.genoss.com

Customer Service & Inquiries

This product includes a 6-month limited warranty. Repairs after this period, or damage caused by user misuse or negligence (even within the warranty period), will be subject to a fee.

Tel. +82-31-888-5100 | Fax. +82-31-8019-9123